

# BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

## Adrenal Protocol

**Indications :** evaluate for adrenal tumor

Scan Type	Helical
Scan Start\ End Location	Above diaphragm to iliac crest
Scan Details	Non contrast (5mm X 5mm) venous phase (5mm X 5mm) 10 minute delay (2mm X 2mm) in abdomen window
IV Contrast	70ml injected 2cc/sec
Scan Delay	65 second venous delay 10 minute delay
Reformat Details	Coronals on venous phase and 10 minute delay
Comments	Pt always drinks oral contrast for this exam, unless specified by radiologist.
Images to send to PACS	Scouts, all 3 phases, and coronals