

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Lumbar Spine

Indications :

Non Contrast: fracture, trauma, disc rupture, disc herniation, stenosis

Scan Type	Helical
Scan Start\ End Location	T12-S1
Scan Details	2mm X 2mm in bone and soft tissue windows
IV Contrast	N/A
Scan Delay	N/A
Reformat Details	3mm X 3mm coronal and sagittal recons in bone window
Comments	Always do angled reconstructions through disc spaces 1-5
Images to send to PACS	Scouts, axials in bone and soft tissue windows, coronals, sagittals, angled axials through the spaces