

# BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

## Adrenal Protocol

**Indications :** evaluate for adrenal tumor

Scan Type	Helical
Scan Start\ End Location	Above diaphragm to iliac crest
Scan Details	Non contrast (5mm X 5mm) venous phase (5mm X 5mm) delay phase (5mm X 5mm)
IV Contrast	70ml injected 2cc/sec
Scan Delay	65 second venous delay 15 minute delay phase
Reformat Details	Recon to 2.5mm on all phases
Comments	Pt always drinks oral contrast for this exam, unless specified by radiologist.
Images to send to PACS	Scouts, all 3 phases, coronals, and recons